

# Wixson Hip Positioner

*Provides stable positioning of a patient during hip surgery*

Designed by  
R.L. Wixson, MD

Product No: 4050

**Optional & Replacement Parts:**

4150-C	[2" Spacer]
4150-C4	[4" Spacer]
4150-EK	[4" Screw for use with 2" Spacer]
4150-EK4	[6" Long Screw for use with two 2" Spacers or one 4" Spacer]
4150-EXT	[2" Spacer with Screw]
4150-EXT4	[4" Spacer with Screw]
4150-9	[9" Post]
4150-12	[12" Post]
4150-PD3	[Set of 3 Small Pads]
4050-LPD	[Large Pad]
4050-BP20	[20" Wide Baseplate]
4050-BP24	[24" Wide Baseplate]

*Video Now*



*Available!*

The Wixson Hip Positioner is used for stable positioning of a patient during total hip and revision surgery. It is designed to be placed on top of the operating table.

The base plate is rubber-backed to reduce slipping on the table. The uprights can easily be slid in and out of the multiple slots in the plate for desired positioning and locked into position with the locking bolt. The complete upright assembly is radiolucent.

The upright pads and the base plate pad are made of semi-dense foam to help prevent pressure points and are sealed with a washable coating. The coating also helps to lessen the possibility of skin breakdown.

The hip positioner consists of: One 10" post with double pads, one 6" post with a single pad, one 20" base plate, one base plate pad, two 2" spacers and one 6" screw for use with larger patients. The pad assembly can be adjusted for additional height and width. The upright posts are modular. The complete unit is radiolucent and autoclavable except for the foam pads and base plate.



**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION




©2007 Innomed, Inc., All Rights Reserved


**Phone** 800 548 2362  
912 236 0000  
**Fax** 912 236 7766  
**Internet** www.innomed.net  
**eMail** info@innomed.net

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404







The upright assembly can be rotated and locked in place.




The upright assembly adjusts by a sliding track to accommodate various sized patients. It is locked in the sliding track by tightening one or two locking bolts.



Spacer assembly.



The pads can be adjusted for height and width.



Stulberg Hip Positioner attaches to table via table clamps.

