

## CLEANING AND STERILIZATION PROCEDURE

# Carbon Fiber Instruments and Foot Piece

### **Recommended hand cleaning instructions:**

- Soak the product in an enzyme solution. Follow the manufacturer's direction for dilution and soaking time.
- Rinse the product in warm tap water.
- Soak or wipe the product down with a hospital approved and/or EPA approved germicide.
- Rinse the product in warm tap water.
- Wash the product with an instrument detergent up to a pH of 9.0
- Rinse the product in warm tap water.
- Dry thoroughly and wrap.

### **Recommended washer sterilizer instructions:**

- Soak the product in an enzyme solution. Follow the manufacturer's direction for dilution and soaking time.
- Put thru washer/sterilizer cycle, continue as normal procedure.

*Please use the following guidelines when sterilizing this product:*

### **Recommended sterilizer instructions:**

- Insure that all parts are thoroughly cleaned and that all movable parts are loose and can move freely.
- Do not place in a milk bath or lubricate.

<b>Sterilizer Type:</b>	<b>Pre-vacuum Steam Sterilizer</b>
Temperature:	133-134 degrees C / 270-272 degrees F
<b>Exposure Time:</b>	<b>4 minutes</b>
Drying Time:	20 minutes
Cooling:	At room temperature

Can be pre-vacuum steam sterilized **wrapped or unwrapped** in two disposable/reusable wraps. Run normal vacuum cycle for your institution.

<b>Sterilizer Type:</b>	<b>Flash Sterilizer</b>
Temperature:	133-134 degrees C / 270-272 degrees F
<b>Exposure Time:</b>	<b>10 minutes</b>
Cooling:	At room temperature

Sterilize unwrapped.

Revised 11/22/11

FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

1.800.548.2362