



Flexible Osteotome System

Provides an assortment of osteotome blades for various orthopedic surgery procedures

- ▶ Sharp, flexible blades are well suited for loosening implants from cement or bony ingrowth fixation
- ▶ Various blade widths and profiles allow great flexibility to follow the implant contours
- ▶ Modular handles are made of high impact surgical stainless steel and have a quick-coupling positive locking mechanism for ease of use and quick blade changes
- ▶ Slap hammer threads into the handle and is designed to facilitate blade removal

PRODUCT NO'S:	
S0011-00	[Complete Set]
Individual instruments included in complete set:	
S1002	[Thin Blade] 3" x 8mm
S1003	[Thin Blade] 3" x 10mm
S1004	[Thin Blade] 3" x 12mm
S1005	[Thin Blade] 3" x 20mm
S1006	[Curved Thin Blade] 3" x 12mm
S1007	[Curved Thin Blade] 3" x 20mm
S1008	[Thin Blade] 5" x 10mm
S1009	[Thin Blade] 5" x 8mm
S1020	[Handle with Quick-Coupling End] 6" <i>One Handle Only with this Product Number</i>
S1133	[Radial Blade] 5" x 10mm
S1120	[Radial Blade] 5" x 12mm (not shown)
S1134	[Radial Blade] 5" x 14mm
S1121	[Radial Blade] 5" x 16mm
S1122	[Radial Blade] 5" x 20mm (not shown)
S2007	[Small Slap Hammer] 12"
9018	[Case] (not shown)

OPTIONAL BLADE PRODUCT NO'S: (Not included in complete set)	
OSTEOTOME BLADES	
S1123	[Extra Long Osteotome Blade] 9" x 8mm
S1135	[Radial Osteotome Blade Medial Curve] 6.75" x 11mm
S1136	[Radial Osteotome Blade Lateral Curve] 6.75" x 11mm
S1137	[Radial Osteotome Blade Medial Curve] 5" x 11mm
S1138	[Radial Osteotome Blade Lateral Curve] 5" x 11mm
CHISEL BLADES	
S1222	[Chisel Blade] 2.5" x 8mm
S1223	[Chisel Blade] 2.5" x 10mm
S1224	[Chisel Blade] 2.5" x 12mm
S1225	[Chisel Blade] 2.5" x 20mm
S1228	[Chisel Blade] 5" x 10mm
S1229	[Chisel Blade] 5" x 8mm
S1230	[Chisel Blade] 5" x 20mm
S1231	[Chisel Blade] 5" x 12mm



FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362

CLEANING AND STERILIZATION PROCEDURE

Flexible Osteotome System

Please use the following guidelines when sterilizing this product:

Sterilizer Type:	Pre-vacuum Steam Sterilizer
Minimum Temperature:	132 degrees C / 270 degrees F
Minimum Cycle Time:	4 minutes
Minimum Dry Time:	10 minutes

Sterilizer Type:	Flash Sterilizer
Minimum Temperature:	132 degrees C / 270 degrees F
Minimum Cycle Time:	4 minutes
Minimum Dry Time:	0 minutes

Innomed recommends that the cleaning and decontamination of instruments follow the guidelines set forth by AORN/HIMA and AAMI. Both physical and chemical (detergent) processes are necessary to minimize the bioburden on all soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris, therefore a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. Carefully inspect hidden areas such as cannulations and recesses to assure any residual materials are removed. Once the items have been cleaned and decontaminated they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization. Innomed recommends the use of a mild enzymatic detergent with a low pH. Do not use multipurpose detergents to wash or soak your instruments. Use a specifically compounded low-suds detergent for all instruments. Detergents designed for surgical instruments, pads and straps are specifically formulated to remove protein, organic debris and blood. The neutral pH balance will not damage stainless steel or tungsten carbide inserts. The solution is gentle enough for manual (hand) as well as ultrasonic cleaning.

INSTRUMENT CARE PROCEDURE

1. Visually inspect instruments before cleaning for cracks, tears and chipped areas.
2. Clean instruments thoroughly after use.
 - a. If you use a pre-soak solution; be certain that it has a neutral pH balance.
 - b. Clean instruments in an open position by either hand or ultrasonic cleaner.
 - c. Use a non-metallic brush (toothbrush) to remove stubborn debris. DO NOT use abrasive cleaning solutions or scouring pads.
 - d. DO NOT expose instruments to bleach.

Detergents designed for surgical instruments are specifically formulated to remove protein, organic debris and blood. The neutral pH balance will not damage stainless steel or tungsten carbide inserts. The solution is gentle enough for manual (hand) as well as ultrasonic cleaning.

Revised 1/1/12

FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362