


Koonin Leg Length Caliper - Small

Designed by Michael Koonin, MD

Designed for use in small incisions to help measure and evaluate pre- and post-THR leg length in conjunction with X-ray calibration and clinical judgement


Utilizes a 5/32" (4 mm) threaded pin (dull on the outer end) in the iliac crest and a 1/8" (3,2 mm) pin in the greater trochanter. (The soft tissue is cleared away and a single drill hole is drilled in the trochanter to accommodate the distal pin, and the hole is marked with methylene blue so it can be easily found.)

Alternatively, a 7.3 mm cannulated screw that accepts a 3.2 mm pin may be used in the greater trochanter. Using the sliding caliper, the difference in leg length measurement before hip dislocation and after the THR procedure helps show the change in leg length.



Hole for a 5/32" (4 mm) threaded pin in the wound just proximal to the acetabulum

Locking screw for caliper



Hole for a 1/8" (3.2 mm) pin in the greater trochanter

Locking screw used on pin at 1st measurement to ensure level with 2nd measurement

PRODUCT NO:

1196

Overall Length: 3.25-4.5" (8,3 cm-11,4 cm)



Scan to Launch Our Website



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362