

Rose Hamstring Tendon Harvester


Designed by Donald J. Rose, M.D., FACS, FAAOS

Designed to easily convert from an open to a closed device without sharp edges to facilitate safe harvesting of hamstring tendon autografts

OPEN

CLOSED

PRODUCT NO:
4692
Overall Length: 15" (38,1 cm)
Shaft Length: 9.5" (24,1 cm)
Internal Diameter 5.5 mm
External Diameter: 8 mm



Harvester placed in open position around isolated hamstring tendon after complete lysis of inferior fibrous bands.

Harvester in closed position capturing tendon, with pes anserinus attachment still intact.

Colinear advancement of harvester, without twisting, separating tendon (under tension) from muscular attachment.

Retrieved tendons. Graft length may be maximized by subsequently avulsing pes anserinus from its tibial attachment by distal traction, after both gracilis and semitendinosus tendons are harvested.

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED [®]



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.