

CLEANING AND STERILIZATION PROCEDURE

Carbon Fiber Instruments and Foot Piece

Recommended hand cleaning instructions:

- Soak the product in an enzyme solution. Follow the manufacturer's direction for dilution and soaking time.
- Rinse the product in warm tap water.
- Soak or wipe the product down with a hospital approved and/or EPA approved germicide.
- Rinse the product in warm tap water.
- Wash the product with an instrument detergent up to a pH of 9.0
- Rinse the product in warm tap water.
- Dry thoroughly and wrap.

Recommended washer sterilizer instructions:

- Soak the product in an enzyme solution. Follow the manufacturer's direction for dilution and soaking time.
- Put thru washer/sterilizer cycle, continue as normal procedure.

Please use the following guidelines when sterilizing this product:

Recommended sterilizer instructions:

- Insure that all parts are thoroughly cleaned and that all movable parts are loose and can move freely.
- Do not place in a milk bath or lubricate.

Sterilizer Type:	Pre-vacuum Steam Sterilizer
Temperature:	133-134 degrees C / 270-272 degrees F
Exposure Time:	4 minutes
Drying Time:	20 minutes
Cooling:	At room temperature

Can be pre-vacuum steam sterilized **wrapped or unwrapped** in two disposable/reusable wraps. Run normal vacuum cycle for your institution.

Sterilizer Type:	Flash Sterilizer
Temperature:	133-134 degrees C / 270-272 degrees F
Exposure Time:	10 minutes
Cooling:	At room temperature

Sterilize unwrapped.

Revised 1/10/17

FREE TRIAL ON MOST INSTRUMENTS

© 2017 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362