

Acetabular Retractors with Large Handle

Designed for retraction around the acetabulum



PRODUCT NO:
C1009 [Standard]
Overall Length: 15.75" (40 cm)
Handle Length: 6.125" (15,6 cm)
Blade Depth: 4.5" (11,4 cm)
Blade Width: 1.1" (28 mm)



PRODUCT NO:
C1010 [90°]
Overall Length: 14.25" (36,2 cm)
Handle Length: 6.125" (15,6 cm)
Blade Depth: 5.25" (13,3 cm)
Blade Width: 1.1" (28 mm)



ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362