



Adjustable Knee & Tibial Positioner

Designed by Ashutosh Chaudhari, MD

Adjustable design allows for use in a multitude of procedures around the knee such as tibial nailing, tibial condyle plating, patella fracture fixation, supracondylar fracture plating, supracondylar fracture nailing, and total knee replacement

Steam sterilizable.



| PRODUCT NO'S: | |
|--------------------------------|---|
| 2770-00 [Set] | Includes Positioner, Pad, and Two Short Straps |
| Individual / Replacement Parts | |
| 2770-01 [Positioner] | Overall Length (Folded): 28" (71,1 cm) Overall Length (Flat): 54.75" (139 cm) Maximum Triangle Height: 14" (35,6 cm) Width: 5.5" (14 cm) Thickness (Folded): 1.8" (4,6 cm) Thickness (Flat): .75" (1,9 cm) |
| 2770-P [Silicone Pad] | Dimensions: 12" x 5.5" (30,5 cm X 14 cm) |
| 2590-S [Short Straps] | Pkg of 10 |



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362

© 2018 Innomed, Inc.