



Amstutz Hip Resurfacing Instruments

Designed by Harlan C. Amstutz, MD

*Amstutz Femoral Head-Neck Elevator
Used for acetabular exposure*



PRODUCT NO'S:
3410 [Wide Blade] Overall Length: 12.25" (31,1 cm) Blade Width at Widest: 67 mm
3410-01 [Narrow Blade] Overall Length: 12.25" (31,1 cm) Blade Width at Widest: 50 mm




*Amstutz Femoral Retractor
Used for acetabular exposure*



PRODUCT NO:
6410 Overall Length: 10.75" (27,3 cm) Blade Width: 18 mm



*Amstutz Charnley-type
Acetabular Exposure Pin Set
Used to enhance exposure
in the acetabulum*



PRODUCT NO'S:
1200-00 [Set: Ins/Ext & Two Pins]
1200-0A [Set: Ins/Ext & Two Pins w/Stop]
Also sold Individually:
1200-01 [Inserter/Extractor]
1200-02 [Pin] Overall Length: 4.5" (11,4 cm) Pin Depth: 2" (5,1 cm) Pin Diam.: 3.9 mm
1200-03 [Pin with Stop] Overall Length: 4.5" (11,4 cm) Pin Tip-to-Stop Depth: .75" (1,9 cm) Pin Diameter: 3.2 mm
1200-04 [Deep Pin] Overall Length: 7" (17,8 cm) Pin Depth: 4.5" (11,4 cm) Pin Diameter: 3.9 mm



Cement Packer & Trimmer



PRODUCT NO:
4995 Overall Length: 9.75" (24,8 cm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362