


# Andrews Modified Tibial Wedge Clamp

Designed by Scott Andrews, MD and Kuldeep Sidhu, MD

*Designed to help remove the cut tibial bone during  
total knee procedures*

The bone is held by the spikes which helps it to come out in one piece, and also helps with release of soft tissues from the bone.



**PRODUCT NO:**

3642  
Overall Length: 10.25"  
Pads: 60mm x 30mm  
Front Spike Length: 14mm  
Back Spike Length: 7.5mm

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362