

Anterior Femoral Punches

Designed by Brandon Thompson, CST/CFA


Designed with a delrin pad to help protect the femoral stem trunnion while removing the femoral head during anterior approach total hip revision arthroplasty



- ▶ Three stem angles allow choice of optimal approach
- ▶ Angled punches allow for better striking force to help break the taper of the head and stem
- ▶ The delrin pad helps prevent scratching of the femoral stem trunnion

THE DELRIN PAD SHOULD NOT BE USED FOR IMPACTION.

PRODUCT NO'S:	
8626-A [Angled Up]	Overall Length: 8.75" (22,2 cm) Up Angle: 40°
8626-L [Left]	Overall Length: 9" (22,9 cm) Left Angle: 40°
8626-R [Right]	Overall Length: 9" (22,9 cm) Right Angle: 40°



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

