


Extension Set for Anterior THR Tables

Designed by David Ott, MD

Designed to add lift to the femoral hook during an anterior THR case and be able to remove without breaking the sterile field



PRODUCT NO'S:
8004-00 [Set of One Each]
Also available individually:
8004-S [Short Extension] Extension Length: 2" (5,1 cm) Overall Length: 2.6" (6,6 cm)
8004-L [Long Extension] Extension Length: 3" (7,7 cm) Overall Length: 3.625" (9,2 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362