


Anterior Watson Jones Total Hip Arthroplasty System

Instrument system specifically designed for Direct Anterior approach THR




PRODUCT NO'S:	
6300-00	[Complete Set]
Included in Set or available Individually:	
6301-L	[Anterior Watson Jones Awl – Left]
1	Overall Length: 12.375" (31,5 cm)
6301-R	[Anterior Watson Jones Awl – Right]
2	Overall Length: 12.375" (31,5 cm)
6302-01	[Anterior Watson Jones Lighted Mueller Retractor]
3	Blade Width: 27 mm Handle Length: 5.5" (14 cm) Overall Length: 12.75" (32,4 cm)
6303-01	[Anterior Watson Jones Wide Lighted Retractor]
4	Blade Width: 40 mm Blade Tip Length: 17 mm Overall Length: 12.75" (32,4 cm)
6304-01	[Anterior Watson Jones Narrow Lighted Retractor]
5	Blade Width: 26 mm Blade Tip Length: 23 mm Overall Length: 16.5" (41,9 cm)
6305	[Anterior Watson Jones 90° Cobra Retractor]
6	Blade at Widest: 37 mm Blade Depth: 7.5" (19,1 cm) Blade Prong: 18 mm Long X 10 mm Wide Overall Length: 12.5" (31,8 cm) Handle Length: 9.5" (24,1 cm)
6306	[Anterior Watson Jones Deep Hohmann Retractor]
7	Blade Width: 17 mm Prong Length: 34 mm Overall Length: 9.25" (23,5 cm)
6307	[Anterior Watson Jones Straight Hohmann Retractor]
8	Prong Length: 4,2 cm Overall Length: 9.5" (24,1 cm)
6308	[Anterior Watson Jones Femoral Starter Drill]
9	Drill Tip Dimensions: 18 mm Long X 12 mm Diameter Overall Length: 6" (15,2 cm) Shaft Length: 5.25" (13,3 cm)

Lighted retractors can be attached to a fiber optic light cable with ACMI (female) connector.




Optional quick-connect driver (NOT INCLUDED IN RETRACTOR SET) for use with the starter drill:

PRODUCT NO:	
8248	[Fixed Driver with Zimmer Hall Quick-connect]
	Overall Length: 5.75" (15,6 cm) Handle Width: 4.625" (11,6 cm)



Fixed Driver with Zimmer Hall Quick-connect



1 Awl – Left #6301-L

2 Awl – Right #6301-R

3 Lighted Mueller Retractor #6302-01

4 Wide Lighted Retractor #6303-01

5 Narrow Lighted Retractor #6304-01

6 90° Cobra Retractor #6305

7 Deep Hohmann Retractor #6306

8 Straight Hohmann Retractor #6307

9 Femoral Starter Drill with Zimmer Hall Quick-connect #6308



Scan to Launch Our Website



FREE TRIAL ON MOST INSTRUMENTS

ISO 13485:2016

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2019 Innomed, Inc.