

Atlatl Super Slap Hammer

Designed for when extra powerful slap hammer force is needed

Repositionable silicone grip handles are available for use with the long version of the Atlatl, and are removable for sterilization.




Slap hammer rod not included - available separately.

PRODUCT NO'S:

3924-L-00 [Long with Silicone Handles] Overall Length: 22" (55,8 cm) Includes (2) 3924-RH Silicone Grip Handles. Slap hammer rod not included.
3924-L [Long without Handles] Overall Length: 22" (55,8 cm) Slap hammer rod not included.
3924-S [Short] Overall Length: 16" (40,7 cm) Slap hammer rod not included.
3924-RH [Silicone Grip Handle] Overall Length: 4" (10,2 cm) (One handle only with this product number)
3925-A [16" Slap Hammer Rod only] (Rod only with this product number)



With the long Super Slap Hammer, two handed force can be used to power the drive weight from bottom to top along the entire length of the slap hammer rod



For use with a 3/8" diameter slap hammer rod, including the Innomed #3925 & #3935 slap hammers on the following extraction instruments:

Hip - Femoral Component

- 3610 Universal Modular Hip Component Extractor - Standard
- 3610-R Universal Modular Hip Component Extractor - Anterior
- 3611 Heck Anterior Modular Hip Component Extractor
- 4175-00 Whelan Hip Stem Extractor
- S1202 Femoral Extraction Instrument - Loop
- S1203 Femoral Extraction Instrument - J-Hook
- S1203 Femoral Extraction Instrument - One-Piece

Hip - Acetabular Cup/Shell/Liner

- 3638-00 Lombardi Hip Cup Liner/Shell Extractor
- 3660 Gorski Hip Cup Extraction Hook - 6.5 mm
- 3665 Gorski Hip Cup Extraction Hook - 5.0 mm

Knee


- 3630 Tibial Knee Component Extractor
- 3920 Femoral Knee Component Extractor
- 3650 4 mm Tibia Tray Removal Hook
- 3655 8 mm Tibia Tray Removal Hook

Shoulder

- 3670 Nicholson Universal Humeral Prosthesis Extractor

General

- 3966 Large Bent Jaw OrthoVise



Hip - Femoral Component

- 3610 Universal Modular Hip Component Extractor - Standard
- 3610-R Universal Modular Hip Component Extractor - Anterior
- 3611 Heck Anterior Modular Hip Component Extractor
- 4175-00 Whelan Hip Stem Extractor
- S1202 Femoral Extraction Instrument - Loop
- S1203 Femoral Extraction Instrument - J-Hook
- S1203 Femoral Extraction Instrument - One-Piece

Hip - Acetabular Cup/Shell/Liner

- 3638-00 Lombardi Hip Cup Liner/Shell Extractor
- 3660 Gorski Hip Cup Extraction Hook - 6.5 mm
- 3665 Gorski Hip Cup Extraction Hook - 5.0 mm

Knee

- 3630 Tibial Knee Component Extractor
- 3920 Femoral Knee Component Extractor
- 3650 4 mm Tibia Tray Removal Hook
- 3655 8 mm Tibia Tray Removal Hook

Shoulder

- 3670 Nicholson Universal Humeral Prosthesis Extractor

General

- 3966 Large Bent Jaw OrthoVise

ISO 13485:2016

Scan to Launch Our Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362

© 2019 Innomed, Inc.