



# Axillary Nerve Protector

Designed by Brett Sanders, MD

*Designed for inferior capsular release during shoulder arthroplasty and glenoid exposure*

The tapered freer end helps separate the axillary nerve and inferior capsule, even in difficult exposures. Non-conductive material allows the use of a bovie knife directly in the small channel cutting guide (on both sides). Reversible for right and left use.



**PRODUCT NO:**

8029

Overall Length: 7.125" (18,1 cm)

Width: 12 mm

Thickness: 4 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

© 2023 Innomed, Inc.