

Bacastow Axillary Nerve Retractor with Suction

Designed by David Bacastow, MD

Designed with a curved tip to slip all the way under the capsule during shoulder surgery, helping to protect the axillary nerve, while also providing suction of smoke away from the surgical site

Made of autoclavable Radel material, the unit is non-conductive of current and resists the high temperatures associated with the use of electrocautery.




PRODUCT NO:

8739

Overall Length: 11" (27,9 cm)

Width: .5" (12,7 mm)

Tongue Length: 1.2" (30 mm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED The Innomed logo features the word "INNOMED" in a bold, sans-serif font, followed by a stylized lowercase "i" that incorporates a human figure. To the right of the logo is a small American flag.

1.800.548.2362