

Bacastow Shoulder Capsular Retractor

Designed by David Bacastow, MD


Designed to help place tension on the inferior capsule for improved visualization and dissection when performing anatomic or reverse shoulder replacement

Rotating arms allow left or right use. Can be used in conjunction with Axillary Nerve Retractor with Suction (#8739) for additional safety.

PRODUCT NO:

5185

Overall Length: 8" (20,3 cm)
Length to Pivot: 5.75" (14,6 cm)
Small Blade Depth: 1" (2,54 cm)
Small Blade Width: 10 mm
Wide Blade Depth: 2" (5,1 cm)
Wide Blade Width: 2.25" (5,7 cm) Maximum
Wide Blade Offset: .75" (1,9 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362