




# Beicker Curette Suction Device

Designed by Clint Beicker, MD

*Designed to help visualization of a fracture site within a fracture hematoma*

Also useful for arthroscopic curettage of osteochondral lesions.



Ultra hard titanium nitride coating helps to extend life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

**PRODUCT NO:**

4231

Overall Length: 10.5" (26,7 cm)  
Curette Cup: 7.5 mm x 5.5 mm



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**