


Bhargava Modified Meniscal Clamp

Designed by Tarun Bhargava, MD

Low-profile design helps facilitate grasping the posterior portion of the meniscus

Improved bite when tension is placed on the meniscus.

Can also be used to help remove the fat pad and suprapatellar bursa.




PRODUCT NO:

1886

Overall Length: 7" (17,8 cm)

Jaw Length: 1.125" (2,9 cm)

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362