


Blair Acetabular Cup Positioner

Designed by Christopher Blair, DO

Designed to help adjust the position of an acetabular cup




PRODUCT NO:

4159

Overall Length: 11.5" (29,2 cm)

Shaft Offset: 1" (2,54 cm)

Head Diameter: 1.18" Inside (30 mm / 1.5" Outside (38 mm))



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362