


# Depth Gauge

*Designed to help provide measurement of the depth/length of any bone hole for proper screw length determination*

**PRODUCT NO:**

8015  
Overall Length: 7.625" (19.4cm)  
Scale: From 0 to 48mm



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71



**1.800.548.2362**