

Bone Graft Impactors

Tap bone graft or bone parts into place with minimal bone trauma

Designed with serrated, stainless steel tips and available in three shapes: round, square and rectangular.



PRODUCT NO'S:	
5310 [Round]	Head Diameter: 12.5 mm Overall Length: 9.5" (24.1 cm) Handle Length: 4.25" (10.5 cm)
5320 [Square]	Head Dimensions: 10 mm x 10 mm Overall Length: 9.5" (24.1 cm) Handle Length: 4.25" (10.5 cm)
5325 [Square with Delrin Tip]	Head Dimensions: 10 mm x 10 mm Overall Length: 9.5" (24.1 cm) Handle Length: 4.25" (10.5 cm)
5330 [Rectangular]	Head Dimensions: 10 mm x 3 mm Overall Length: 9.5" (24.1 cm) Handle Length: 4.25" (10.5 cm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362