


Boynton Punch

Designed by L. Boynton, MD

*Helpful in removing trial, femoral
and revision total knee components*


The flange end fits onto the flange of a femoral
knee component or trial.



PRODUCT NO'S:

5120-01 [Standard]
Overall Length: 11.75" (29,8 cm)
Shaft Diameter: 9,5 mm

5120-02 [Offset]
Overall Length: 11.75" (29,8 cm)
Shaft Diameter: 9,5 mm
Punch End Offset: 60 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.