

# Bozeman Cement Trimmer

Designed by Daniel M. Gannon, MD

*Combines the two most common cement trimming tools combined into one*

The tool has a blunt blade tip on one end to help with separation of the trimmed cement. The angled curette end helps gather the trimmed cement. The thin shank and angled curette can reach into tight spaces such as the back of the implants to remove excess cement. The ends are titanium nitrite coated to help eliminate metal transfer.



**PRODUCT NO:**  
5245  
Overall Length: 8.5" (21,6 cm)

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

© 2023 Innomed, Inc.