

Bradley Periosteal Elevator

Designed by Gary W. Bradley, MD

Ultra hard titanium nitride coating helps to extend blade life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.



PRODUCT NO'S:

4719 [1/2"]

Overall Length: 11" (27,9 cm)

Blade Width: .5" (13 mm)

4720 [3/4"]

Overall Length: 11" (27,9 cm)

Blade Width: .75" (19 mm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362