




Browne Deltoid Retractors

Used for the Delto-Pectoral approach, the retractor contours the humeral head with effortless deltoid retraction allowing extensive exposure



Small

Small Deep

Large



The small deep version is designed for use in large patients.

PRODUCT NO'S:	
1670-01 [Small]	Blade Width: 45 mm Overall Length: 11.5" (29,2 cm)
1670-01D [Small Deep]	Blade Width: 45 mm Overall Length: 12.25" (31,1 cm)
1670-02 [Large]	Blade Width: 57 mm Overall Length: 11.5" (29,2 cm)

MADE FOR INNOMED IN GERMANY

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362