


# Budny Wire Drill Guide

*Designed to provide stability for the insertion of smooth and olive wires during the application of ring-based external fixation systems*


Designed by  
Adam M. Budny, DPM

Product No: 1188  
*Wire not included.*



The 10cm long cannula helps to provide stability and directional control during wire insertion

Pistol grip allows ambidextrous use



Internal stainless steel cannula (2mm) can be completely exposed by opening the swinging latch-cover mechanism, allowing easy insertion/release of a wire

The entire unit is autoclavable.

**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**  
912 236 0000  
Fax 912 236 7766  
Internet [www.innomed.net](http://www.innomed.net)  
eMail [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

