

Byrd Arthroscopic Hip Elevators

Designed by J. W. Thomas Byrd, MD

Designed for use in various arthroscopic hip procedures

Versatile **elevators** allow freeing up of soft tissue around bone, can assist in removing the fibrous layer off of a cam lesion, and can be used to help elevate the capsule. They are also useful for mobilizing the labrum in preparation for repair or acetabular rim trimming.



PRODUCT NO'S:
Overall Length: 13.625"
2045 [Extra Long Elevator – Square Tip] Pad Dimensions: 4.8mm x 19mm
2046 [Extra Long Elevator – Angled, Joseph Straight Edge] Pad Dimensions: 4.4mm x 12mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362