



Calcaneal Spreader

Designed by Michael Forness, DO

Separates the calcaneal osteotomized bone for placement of tricortical bone graft

Pads have a large surface area, which easily separates the calcaneal osteo-tomized bone for placement of tricortical bone graft. Large pad surface area helps prevent the compression of soft calcaneal cancellous bone.



Smooth Pads



Grooved Pads

PRODUCT NO'S:
1880 [Smooth Pads] Overall Length: 7" (17,8cm) Pad Dimensions: 15mm x 12mm
1881 [Grooved Pads] Overall Length: 7" (17,8cm) Pad Dimensions: 15mm x 12mm

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362