

Calibrated Ortho Spreader with Slotted Tips

Designed by Jason Bariteau, MD

*A lamina spreader with a very thin closed profile,
designed to enable distraction in tight joint spaces like
the subtalar and talonavicular joints*



PRODUCT NO:

1841

Overall Length: 6.75" (17.1 cm)

Prong Length: .5" (12.7 mm)

Calibrations: 10 mm to 35 mm



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362