

# Calibrated Ortho Spreader with Slotted Tips

Designed by Jason Bariteau, MD

*A lamina spreader with a very thin closed profile,  
designed to enable distraction in tight joint spaces like  
the subtalar and talonavicular joints*



**PRODUCT NO:**

1841

Overall Length: 6.75" (17.1 cm)

Prong Length: .5" (12.7 mm)

Calibrations: 10 mm to 35 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

# INNOMED