

Calibrated Ortho Spreader without Teeth

Helps separate the femur and tibia during knee replacement procedures

Available with flat or serrated outside blades. Calibrated ratchet (in mm) is helpful in ligament balancing.

- ▶ Can also be used for spine surgery where the calibrated ratchet (in mm) is used to accurately measure the size of opening - useful in procedures to help assess bone graft needs.
- ▶ Also used for foot & ankle surgery

PRODUCT NO'S:	
Flat Outside Pads	Serrated Outside Pads
1842 [Small Flat] Overall Length: 6.5" Blade Width: 7mm Blade Thickness: 1.68mm	1842-01 [Small Serrated] Overall Length: 6.5" Blade Width: 7mm Blade Thickness: 1.68mm
1843 [Medium Flat] Overall Length: 9.25" Blade Width: 10mm Blade Thickness: 1.68mm	1843-01 [Medium Serrated] Overall Length: 9.25" Blade Width: 10mm Blade Thickness: 1.68mm

**PROUDLY
MADE
IN THE
USA**



Calibrated ratchet (in mm) is helpful in ligament balancing.

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362