

Calibrated Ortho Spreader without Teeth

Calibrated ratchet (in mm) is used to accurately measure the size of opening - useful in procedures to help assess bone graft needs.

- ▶ Can also be used to help separate the femur and tibia during knee replacement procedures
- ▶ Also used for foot & ankle surgery

PRODUCT NO'S:	
Flat Outside Pads	Serrated Outside Pads
1842 [Small Flat] Overall Length: 6.5" Blade Width: 7mm Blade Thickness: 1.68mm	1842-01 [Small Serrated] Overall Length: 6.5" Blade Width: 7mm Blade Thickness: 1.68mm
1843 [Medium Flat] Overall Length: 9.25" Blade Width: 10mm Blade Thickness: 1.68mm	1843-01 [Medium Serrated] Overall Length: 9.25" Blade Width: 10mm Blade Thickness: 1.68mm

PROUDLY
MADE IN
THE USA



Available with
flat or serrated
outside blades.

Calibrated ratchet (in mm) is used to accurately measure the size of opening - useful in procedures to help assess bone graft needs.

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362