



# Cannestra Trochanteric Fracture Reduction Clamp

*Designed to help reduce comminuted intertrochanteric and subtrochanteric hip fractures, this clamp is offset at its ends to avoid placement into the fracture bed*



Clamping ends are curved and rotated to allow maximum bony contact upon fracture reduction. Ideal for fractures with a flexed anterior cortical spike. Made for right and left hip fracture configurations.

Product No's:  
3860-L [Left]  
3860-R [Right]

OVERALL LENGTH: 11.25"

**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**  
912 236 0000  
Fax 912 236 7766  
Internet [www.innomed.net](http://www.innomed.net)  
eMail [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

