

Cartilage Graspers

Helps to grasp and hold cartilage, tendons, soft tissues and loose bodies

Shaft allows for use in narrow spaces.



Saw Tooth Jaw



Saw Tooth Jaw with 6" (15,2 cm) Shaft.


PRODUCT NO:

1785 [Saw Teeth]
Shaft Length: 6" (15,2 cm)
Overall Length: 9.25" (23,5 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



Shark Tooth Jaw



Shark Tooth Jaw with 8" (20,3 cm) Shaft.

Shark Tooth Jaw with 5" (12,7 cm) Shaft.

Shark tooth designed by Luis Ulloa,
design modification by Michael Soudry, MD

PRODUCT NO'S:

1777 [5" with Shark Teeth]
Shaft Length: 5" (12,7 cm)
Overall Length: 8.25" (21 cm)
Jaw Bite: 2 mm x 6.5 mm

1779 [8" with Shark Teeth]
Shaft Length: 8" (20,3 cm)
Overall Length: 11.25" (28,6 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362