



# Chandler Retractor

*Used for retracting tissue away from the bone, and helpful for posterior exposure of the tibia in MIS surgery*

Allows the surgeon to retract soft tissue away from the bone, and can be used for hip and knee surgery. The handle is contoured away from the field of view and working area. Available in three blade sizes: 5/8", 3/4" and 1".

The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.



PRODUCT NO'S:	
3220-01 [5/8"]	Overall Length: 9.875" (25,1 cm) Blade Width: 5/8" (1,6 cm)
3220-02 [3/4"]	Overall Length: 9.875" (25,1 cm) Blade Width: 3/4" (1,9 cm)
3220-04 [1"]	Overall Length: 9.875" (25,1 cm) Blade Width: 1" (2,54 cm)
3220-02R* [OrthoLucent™ 3/4"]	Overall Length: 9.875" (25,1 cm) Blade Width: 3/4" (1,9 cm)



\* MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND

[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.