

Chandran Anterior Retractor for THR

Designed by Rama E. Chandran, MD

Design helps to expose the anterior rim of the acetabulum and helps prevent displacement of the retractor while reaming the acetabulum during direct anterior hip replacement



PRODUCT NO:
6311
Overall Length: 13.5" (34,3 cm)
Depth from Flat: 4.75" (12,1 cm)
Blade Width: .625" (1,6 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED