



# Chandran Bent Serrated Curette

Designed by Rama E. Chandran, MD

*Serrated design allows for easier removal of cancellous bone in the proximal femur in total joint arthroplasty*

**PRODUCT NO:**

5171  
Overall Length: 11.75" (29.8 cm)  
Handle Length: 5.5" (14 cm)  
Cup Size: 7 mm X 12 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71



**1.800.548.2362**