

# Chandran Distal Biceps Tissue Protector

Designed by Rama E. Chandran, MD

*Designed to help protect tissue and expose the radial tuberosity during distal biceps tendon repair*

Using downward pressure, the teeth help to engage bone to keep the protector in place.

*Also useful to help expose the humerus during proximal subpectoral biceps repair.*

**PRODUCT NO:**


3224

Overall Length: 8" (20,3 cm)

Handle Length: 5.5" (14 cm)

Blade Depth: 1.7" (4,3 cm)

Blade Width: .65" (1,7 cm)



*Helps to protect tissue and expose the radial tuberosity during distal biceps tendon repair.*

*Also useful to help expose the humerus during proximal subpectoral biceps repair.*



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**