



Chandran Modified Knee Retractor

Designed by Rama E. Chandran, MD

Teeth designed to help prevent tilting of the retractor and protect the patellar tendon during robotic assisted total knee replacment

Also useful to retract structures on the lateral side of the tibia.



PRODUCT NO:
7117
Overall Length: 8.75" (22,2 cm)
Blade Width: .625" (1,6 cm)
Blade Depth from Handle: 5" (12,7 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362