


Chandran Tibial Knee Retractor

Designed by Rama E. Chandran, MD

Designed for use in TKR, the hook on the front of the blade acts as a stop to help prevent the retractor from deep penetration behind the tibia

PRODUCT NO:

4533
Overall Length: 10.5" (26,7 cm)
Depth from Bend: 3" (7,6 cm)



ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362