

# Charnley Type Tissue Needle Forceps

Designed by Amal Das Jr., MD

*Helpful for wound closure in deep areas with fascia under tension, the forceps can help retrieve a needle in a tight spot*

**PRODUCT NO:**

1165

Overall Length: 6.875"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



INNOMED 

1.800.548.2362