


Chung T-Handle Retractors

Designed by Raymond Chung, MD


Designed with a T-handle for easier holding and to help reduce finger and thumb fatigue

PRODUCT NO'S:


1159 [Sharp Rake] Overall Length: 4.5" (11,4 cm) Blade Width: 9 mm Blade Depth: 7 mm	1159-01 [Extended Sharp Rake] Overall Length: 5.625" (14,4 cm) Blade Width: 9 mm Blade Depth: 7 mm	New!
1161 [Blunt Rake] Overall Length: 4.5" (11,4 cm) Blade Width: 9 mm Blade Depth: 7 mm	1161-01 [Extended Blunt Rake] Overall Length: 5.625" (14,4 cm) Blade Width: 9 mm Blade Depth: 7 mm	New!
1162 [Senn] Overall Length: 4.5" (11,4 cm) Blade Width: 6 mm Blade Depth: 16 mm	1162-01 [Extended Senn] Overall Length: 5.625" (14,4 cm) Blade Width: 6 mm Blade Depth: 16 mm	New!



Standard Shaft



Extended Shaft



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362