



Coated Allis Bone Clamps

Modification of design by Charles T. Resnick MD

A traditional Allis Bone Clamp designed with a longer ratchet—for a wider opening to allow a bone and plate to be clamped and locked onto—and coated end(s) to prevent from marring a component surface



PRODUCT NO'S:	
1381	[One Coated End]
Overall Length: 6.125" (15,9 cm)	
Ratcheted Clamp Opens to: 35 mm	
Non-coated-end Width: 4 mm	
1382	[Two Coated Ends]
Overall Length: 6.125" (15,9 cm)	
Ratcheted Clamp Opens to: 35 mm	
Non-coated-end Width: 4 mm	



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2017 Innomed, Inc.