



Cobra Retractors

A general purpose instrument for use around the femur and acetabulum

Standard with Sharp Tip



Standard



Standard with Ergonomic Handle



Lighted Standard



OrthoLucent™ Standard



Medium



Wide



Lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector. Retractor can be steam sterilized.

The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

PRODUCT NO'S:

6129 [Standard w/Sharp Tip] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
6130 [Standard] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
6130-H [Standard with Ergo. Handle] Overall Length: 12" (30,5 cm) Ergonomic Handle Length: 5" (12,7 cm) Blade at Widest: 33 mm
6130-L [Lighted Standard] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6130-R* [OrthoLucent™ Standard] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6132 [Medium] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 46 mm
6140 [Wide] Overall Length: 11.75" (29,8 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 56 mm



* MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND



Light Source Cable Adapters

PRODUCT NO'S:

8009-S [ACMI to Storz Adapter]
8009-W [ACMI to Wolf Adapter]



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2018 Innomed, Inc.