



Cobra Retractors

A general purpose instrument for use around the femur and acetabulum

Standard with Sharp Tip
#6129



Standard
#6130



Lighted Standard
#6130-L



OrthoLucent™ Standard
#6130-R



Standard with Blunt Teeth
#6130-01



Round and square tip with additional teeth on the ends for better gripping



Straight Tip with Blunt Teeth
#6130-02



Medium
#6132



Wide
#6140



The lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector and can be steam sterilized.

PRODUCT NO'S:

6129 [Standard with Sharp Tip] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
6130 [Standard] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
6130-L [Lighted Standard] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6130-R* [OrthoLucent™ Standard] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6130-01 [Standard with Blunt Teeth] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
6130-02 [Straight with Blunt Teeth] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
6132 [Medium] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 46 mm
6140 [Wide] Overall Length: 11,75" (29,8 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 56 mm

New!



The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.



* MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.