



# Cobra Retractors

*A general purpose instrument for use around the femur and acetabulum*

**Standard with Sharp Tip**  
#6129



**Standard**  
#6130



**Lighted Standard**  
#6130-L



**OrthoLucent™ Standard**  
#6130-R



**Standard with Blunt Teeth**  
#6130-01

*New!*



*Teeth help to prevent slippage*

**Straight Tip with Blunt Teeth**  
#6130-02

*New!*



**Medium**  
#6132



**Wide**  
#6140



The lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector and can be steam sterilized.

The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

## PRODUCT NO'S:

6129 [Standard with Sharp Tip]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade at Widest: 33 mm

6130 [Standard]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade at Widest: 33 mm

6130-L [Lighted Standard]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade Width at Widest: 33 mm

6130-R\* [OrthoLucent™ Standard]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade Width at Widest: 33 mm

6130-01 [Standard with Blunt Teeth]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade at Widest: 33 mm

*New!*

6130-02 [Straight with Blunt Teeth]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade at Widest: 33 mm

*New!*

6132 [Medium]  
Overall Length: 12" (30,5 cm)  
Handle Length: 7" (17,8 cm)  
Blade at Widest: 46 mm

6140 [Wide]  
Overall Length: 11.75" (29,8 cm)  
Handle Length: 7" (17,8 cm)  
Blade at Widest: 56 mm



\* MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND

[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED**



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.