



# Collateral Ligament Retractor

*Helps protect the lateral collateral ligament while exposing the proximal tibia*


Used during total knee surgery and is inserted between the lateral collateral ligament and bone to protect the ligament and expose the proximal tibia. The dual prongs keep the retractor from rocking and assist in the insertion. The retractor is bent so that it is out of the way of the operating surgeon.

**PRODUCT NO:**

6620

Overall Length: 8" (20,3 cm)

Prong Width: 5 mm | 11 mm Gap | 5 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

© 2023 Innomed, Inc.