



Craig-Type Pin Extractor Set

Designed to firmly tighten circumferentially around a wire, pin, broken screw, etc. for removal — especially helpful for the removal of threaded pins

- Removes pins & screws up to 5mm (.2") diameter and wires as small as 1/32" (.03") diameter
- Five interchangeable collets for various grasping capacities
- Two cross-handle insert rods give strong leverage for locking the collet securely onto the pin





Slap Hammer

Handle Draw Bar

Closing Sleeve with Knurled Hand Wheel

Collets (5 Sizes)

Cross-handle Insert Rods



Product No: 1215-00
Includes Sterilization Case

FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000

Fax 912 236 7766

Internet www.innomed.net

eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

