

Curved Canal Rasps

Design modification by Michael Messieh, MD of original design by Anthony Unger, MD.

Designed for preparation of the femoral canal for insertion of a cemented or cementless hip stem, the multiple diameters serve to prepare the femoral canal after the initial 5 mm is used to find the curvature of the canal



PRODUCT NO'S:

3004-01-08 [8 mm]
Overall Length: 11" (27,9 cm)
Handle Length: 5" (12,7 cm)

3004-01-10 [10 mm]
Overall Length: 11" (27,9 cm)
Handle Length: 5" (12,7 cm)

3004-01-12 [12 mm]
Overall Length: 11" (27,9 cm)
Handle Length: 5" (12,7 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362